

# NIST/EPA/NIH EI-MS LIBRARY

## 2023 Release

### 40K NEW NIST MEASURED/EVALUATED COMPOUNDS

#### 394K Electron Ionization (EI) Spectra

- 347,100 Compounds, 46,954 Replicate Spectra
- 40 K More Compounds than NIST 20


#### 492K Retention Index (RI) Values

- 153K RI Compounds with EI, >40 K Increase
- AI-RI Estimates for All EI Compounds


### COMPOUNDS SELECTED FOR ANALYTICAL RELEVANCE

FULLY EVALUATED WITH DERIVATIVES, RETENTION INDICES & CLASS INFORMATION

| Citation               | New  | Total |
|------------------------|------|-------|
| Wikipedia              | 1570 | 6227  |
| EPA Tox                | 2969 | 7117  |
| Food DB                | 582  | 4273  |
| EU Contaminants        | 6263 | 15149 |
| Protein Data Bank      | 1194 | 4716  |
| Human Metabolite DB    | 1992 | 9393  |
| PFAS                   | 161  | 749   |
| Adams (Essential Oils) | 2136 | 2136  |



Plant Stilbenoid



Mycotoxin ENN A

#### Compound Data


**Name:** Valeryl fentanyl  
**Formula:** C<sub>24</sub>H<sub>32</sub>N<sub>2</sub>O  
**MW:** 364 **Exact Mass:** 364.251463 **CAS#:** 122882-90-0 **NIST#:** 463844 **ID#:** 18711  
**Contributor:** NIST Mass Spectrometry Data Center  
**InChIKey:** VCCPXHWAJYWQMR-UHFFFAOYSA-N **Non-stereo**  
**Synonyms:**  
 1 Pentanamide, N-phenyl-N-[1-(2-phenylethyl)-4-piperidinyl]-; 2 N-Phenyl-N-[1-(2-phenylethyl)-4-piperidinyl]-  
**Other DBs:**  
 - Wikipedia  
 - wikipedia2021\_via\_cheminfo\_SMILES  
 - Environmental  
 - SUSDATAFY22  
**Experimental RI median±deviation (#data)**  
 Semi-standard non-polar: 2962±18 (2)  
 Standard non-polar: 2958±N/A (1)  
**Estimated non-polar retention index (n-alkane scale):**  
 Value: 2938 iu  
 Confidence interval (Nitrogen-containing): 83(50%) 356(95%) iu  
**Retention index:**  
 1. Value: 2980.5 iu  
 Column Type: Capillary  
 Column Class: Semi-standard non-polar  
 Active Phase: HP-5MS  
 Column Length: 30 m

#### Collections

#### RI Averages

#### AI-RI Estimates

#### RI Measured



Piperitol with TMS/Me Derivatives (TFA, Acetyl not shown)

### ENHANCEMENTS

### Hit List

### Compound/Spectrum Lookup

**New: Retention Index For All Hits (Expt'l or AI-RI\*)**

**Select Columns Mass Difference, Reverse Score, Probability, ...**

**New: Expanded Compound Classes**

| #  | Lib. | Match | R.M. | RI    | Del... | S | DBs     | Name                           |
|----|------|-------|------|-------|--------|---|---------|--------------------------------|
| 1  | R    | 999   | 999  | 1442  | 0      | 1 |         | 2-Aminophenol, 2TMS deriva     |
| 2  | M    | 965   | 965  | 1930  | -50    | 0 |         | 3-Amino-2-naphthol, N,O-bis    |
| 3  | M    | 965   | 965  | 1489* | -14    | 1 |         | 2-Amino-m-cresol, N,O-bis(tr   |
| 4  | R    | 956   | 956  | 1309  | -1     | 9 | 38 W... | Catechol, 2TMS derivative      |
| 5  | M    | 955   | 957  | 1582  | -17    | 1 | 5 E     | 4-Mercaptophenol, 2TMS deri    |
| 6  | M    | 953   | 953  | 2023* | -49    | 1 |         | 2,3-Naphthalenediamine, 2TI    |
| 7  | M    | 949   | 949  | 1449  | -29    | 3 | 7 EFGN  | 3-Ethylcatechol, 2TMS          |
| 8  | M    | 947   | 947  | 1600  | -34    | 4 |         | 2-amino-5-chlorophenol, N, C   |
| 9  | M    | 941   | 942  | 1388  | -15    | 2 | 21 W... | 4-Methylcatechol, 2TMS deri    |
| 10 | M    | 939   | 951  | 1778  | -33    | 0 | 1 M     | 1,4-Benzenedithiol, S,S'-bis(t |
| 11 | M    | 938   | 946  | 1469  | -17    | 1 | 7 EM    | 2-Mercaptophenol, 2TMS de      |
| 12 | M    | 936   | 945  | 1755  | -33    | 1 | 2 E     | 1,3-Benzenedithiol, 2TMS de    |
| 13 | M    | 933   | 937  | 1470  | -43    | 9 | 8 EGM   | 3-Isopropyl-1,2-benzenediol,   |
| 14 | M    | 931   | 931  | 1550  | -28    | 1 |         | 2-Amino-4-ethylphenol, 2TMS    |
| 15 | M    | 929   | 929  | 1930  | -50    | 0 |         | 2-Amino-1-naphthol, N,O-bisi-  |

#### Names

#### Spectra Replicates Derivatives Stereo

# NIST EI LIBRARY SOFTWARE


## 2023 EI Release EI MS ANALYSIS TOOLS

### NISTMS

FULL FEATURED MS LIBRARY SEARCH/DISPLAY PROGRAM

MULTIPLE SEARCH TYPES & DISPLAY MODES

5 VIEWS: SPECTRUM SEARCH, FEATURE SEARCH, COMPARE, NAME/SPECTRUM, USER LIBRARY



- 1 Query spectrum list
- 2 Score Histogram
- 3 Hit List –multiple values
- 4 Query spectrum
- 5 Query/Spectrum Compare
- 6 Library Spectrum

5 Tab Views


### Hybrid Search

FOR COMPOUNDS NOT FOUND IN LIBRARY & ID CONFIRMATION

FINDS 'MODIFIED' LIBRARY IDS AND MASSES OF MODIFICATIONS WITH THEIR SHIFTED PEAKS

USES MW ESTIMATE

### DELTA MASS => CHEMICAL FORMULA




| DMass | Replace or Insert          |
|-------|----------------------------|
| 1     | H->D, C->C13, NH->O, CH->N |
| 2     | CH2->O, C=C->C-C           |
| 12    | CH2->C=CH2                 |
| 14    | X-Y->X-CH2-Y               |
| 16    | X-Y->X-O-Y                 |
| 17    | NH->S                      |
| 18    | H->F                       |
| 28    | X-Y->X-CO-Y                |
| 30    | H->CH3O-H                  |
| 32    | X-Y->X-S-Y                 |
| 34    | H->Cl                      |
| 56    | Phenyl->Naphthyl           |
| 70    | H->Phenyl                  |
| 162   | H->Glucose                 |

### MS Interpreter


CONNECT PEAKS TO PLAUSIBLE FRAGMENTS (IN RED)

CONFIRM ID  
COMPUTE FRAGMENT MASSES  
CONNECT PEAKS TO STRUCTURES



### AMDIS

'PURIFIES'  
SPECTRA AND  
CONNECTS  
TO NISTMS



### NIST Website chemdata.nist.gov

#### Libraries, Tools, Service



This site provides information and access to NIST mass spectral data products. A variety of data products are available, including EI and tandem MS libraries (small molecules and peptides), a GC retention index collection as well as certain freely available, specialized spectral libraries. Freely available data analysis tools include AMDIS (Automated Mass Spectral Deconvolution and Identification System for GC/MS), the Mass Spectrum Interpreter (connects chemical structures with mass spectra), and the Mass Spectral Digitizer Program. Also available is a fully functional, version of NIST's Web Search Program v2.39 with a small

#### Tools

- Mass Spectrum Interpreter – Major New Release – February 2019 (v. 3.4). Information and downloads for version 3.4 of this program which connects mass spectral peaks to their probable chemical structure origin (EI and MS/MS, both nominal and accurate mass).
- NIST MS Software and Data - updates, demo, documentation, MSPepSearch, Lib2NIST, RUS libraries and support programs.
- AMDIS – computer program that extracts spectra for individual components in a GC/MS data file (Instructions for using AMDIS with MS Search – 11-25-2019)
- Mass Spectrum Digitizer Program – a tutorial on how to use the program (includes program download) that allows the digitization of graphical spectra
- The NIST Glyco Mass Calculator – a tool to aid in the analysis of glycoforms
- DIMEDR - A Novel Algorithm for Agglomerating Incongruent LC-MS Metabolomics Datasets.
- MS\_Plano (New, 2021) - A new software tool for annotating peaks in collision induced dissociation (CID) tandem mass spectra of peptides and N-glycopeptides.

#### Recurrent Unidentified Spectral Libraries

There are three NIST user libraries of recurrent unidentified spectra (RUS):

- **Food:** A set of 650+ spectra extracted from a set of dried food material, some of these spectra have tentative identifications. These experiments were done with methoximation and TMS derivatization. Data
- **PedUrine:** A set of 200+ spectra from a large set of pediatric urine samples. All of these samples were derivatized with TMS after forming the ethylxime for the non-acid carboxylic groups; the majority of this
- **EssOil:** A set of 1000+ spectra derived from a large set of essential oils (both commercial and laboratory distilled), solvent extract of various plant materials (leaves, flowers, roots, etc). Most of these data were



<http://chemdata.nist.gov>

# NIST Tandem Mass Spectral Library

## 2023 Release

**51,501** Compounds, **60%** More than 2020  
**400 K** Precursor Ions – **2.4 M** Spectra





### Fragmentation Methods

49,590 HRAM (High Res Accurate Mass) Compounds  
50,071 QTOF, HCD, IT-HRAM, QqQ Compounds  
49,561 Ion Trap Compounds (Low Res., up to MS<sup>4</sup>)  
561 APCI HRAM Compounds

### Precursor Ion Types

44,191 Protonated  
19,620 Deprotonated  
14,318 Water/Ammonia Loss  
44,547 Other In-Source Generated

## ALL COMPOUNDS SELECTED FOR RELEVANCE MEASURED AT NIST, THOROUGHLY EVALUATED



| Citation            | New  | Total |
|---------------------|------|-------|
| Wikipedia           | 1618 | 6424  |
| EPA Tox             | 3181 | 8146  |
| Food DB             | 602  | 4491  |
| EU Contaminants     | 6553 | 15818 |
| Protein Data Bank   | 1246 | 4945  |
| Human Metabolite DB | 2071 | 9686  |
| PFAS                | 90   | 116   |

### All Spectra Shown for Each Compound

1. Compound Name  
Many Synonyms

2. All Related Spectra  
All Energies  
Fragmentation Types  
In-Source Ions, ...

Screenshot of the NIST MS Search software interface. The search term is 'Forsyoside'. The results show a list of related compounds and their synonyms. The selected compound is Forsyoside, with a mass spectrum displayed. The spectrum shows a base peak at m/z 267.1 and other significant peaks at m/z 181.1, 203.1, 225.1, 247.1, and 269.1. The chemical structure of Forsyoside is shown above the spectrum.

Full Spectrum  
Selected from List  
at Lower Left

All Spectrum Text  
Information

### Hybrid Search Recommended for All Tandem MS Searches

Find Related  
Compounds Even if  
Query is not in Library

Differences in masses  
are 'Modifications'

Screenshot of the NIST MS Search software interface showing a hybrid search for 'Fentanyl'. The search results show a list of related compounds and their synonyms. The selected compound is Fentanyl, with a mass spectrum displayed. The spectrum shows a base peak at m/z 337.2774 and other significant peaks at m/z 188.1426, 229.2284, 251.2506, 273.2728, 295.2950, 317.3172, 339.3394, 361.3616, 383.3838, 405.4060, 427.4282, 449.4504, 471.4726, 493.4948, 515.5170, 537.5392, 559.5614, 581.5836, 603.6058, 625.6280, 647.6502, 669.6724, 691.6946, 713.7168, 735.7390, 757.7612, 779.7834, 801.8056, 823.8278, 845.8500, 867.8722, 889.8944, 911.9166, 933.9388, 955.9610, 977.9832, 999.0054. The chemical structure of Fentanyl is shown above the spectrum.

Query Spectrum

Shifted Library Peaks  
(gray->red)  
Contain Modification

Library Spectrum

# NIST TANDEM LIBRARY SOFTWARE


## 2023 Tandem Release

### TANDEM (MS/MS) ANALYSIS TOOLS

#### NISTMS

ALL SPECTRA FOR EACH COMPOUND IN ONE LIST

MULTIPLE ENERGIES ION MODES FRAGMENTATION TYPE IN SOURCE IONS



#### FILTER IDS BY SPECTRUM TYPE

M/SMS Hit List Filter

Enable Filtering (Tandem Only)

Instrument Type:  Ion Trap,  Ion Trap + FT,  Other

Allow Precursors:  C16-H,  C16-H + Na, K, Cl, HCO2,  H2O, H4O,  All other (demon/losses),  Exclude isotopic precursors

Preferred HCE/Voltage: 20

Priority: Any


Buttons: OK, Cancel, Help

#### Hybrid Search

-FINDS COMPOUNDS NOT IN LIBRARY AND CONFIRMS IDS  
-USES PEAKS AND LOSSES  
-SHOWS SHIFTED PEAKS

RECOMMENDED FOR ROUTINE USE

| # | L. | S... | D... | DeltaMa... | DBs      | Prec. ... | Inst... | En... | Name                           |
|---|----|------|------|------------|----------|-----------|---------|-------|--------------------------------|
| 1 | hr | 935  | 999  | -0.0000    | 20 EGM   | [M-H]     | HCD     | 20%   | Perfluoroheptanoic acid [M-H]  |
| 2 | hr | 891  | 976  | -49.9968   | 27 WC... | [M-H]     | HCD     | 20%   | Perfluoroheptanoic acid [M-H]  |
| 3 | hr | 866  | 938  | -99.9936   | 24 WC... | [M-H]     | HCD     | 14%   | Perfluorononanoic acid [M-H]   |
| 4 | hr | 511  | 796  | -249.984   | 18 CE... | [M-H]     | HCD     | 27%   | Perfluorododecanoic acid [M-H] |
| 5 | hr | 466  | 714  | -149.990   | 22 WC... | [M-H]     | HCD     | 9%    | Perfluorodecanoic acid [M-H]   |
| 6 | hr | 460  | 828  | -199.987   | 18 CE... | [M-H]     | HCD     | 44%   | Perfluoroundecanoic acid [M-H] |



PFAS ILLUSTRATION: APPLICABLE TO MANY CLASSES

#### MS Interpreter

DIRECT ACCESS FROM NISTMS


CONNECT PEAKS TO STRUCTURE

#### Chemical Formula Calculator

#### Isotope Calculator

Any resolution +/- Charge

Multiple Display Options



Formula Calculator

Mass: 1.031634, mDa: 500

Formula: C-2 H-3 N-1 O-1 S-1

|        |     |     |         |      |   |   |   |   |   |
|--------|-----|-----|---------|------|---|---|---|---|---|
| 1 Ion  | O+E | RDB | Mass    | mDa  | C | H | N | O | S |
| H3NO-1 | Odd | 0   | 1.03109 | -0.5 | 0 | 3 | 1 | 1 | 0 |

Find Formula from Delta Mass

Select Peak Show Structure Fragment

#### NIST Website Chemdata.nist.gov

FREELY AVAILABLE LIBRARIES AND SOFTWARE

#### Downloadable Libraries

Annotated Recurrent Unidentified Spectra  
Urine, Plasma/Serum (ARUS)  
Tryptic Peptides (Human, ...)  
Oligosaccharide Libraries (Milk)  
Glycopeptides (mAb, Glycan Distributions)  
Acyl Carnitines

#### Software

NISTMS (2019)  
MS Piano (Peptide a Annotation)  
MS Pepsearch – General Search Utility for NIST Libraries  
MS Interpreter  
Lib2NIST – Library Conversion



<http://chemdata.nist.gov>